

CONTINUOUS FLOW TRANSFECTION WITH FLOWFECT® TECHNOLOGY ENABLES RAPID OPTIMIZATION AND SCALE-UP FOR EFFICIENT CO-DELIVERY OF mRNA AND PLASMID DNA IN PRIMARY HUMAN T CELLS

Maria Lai, Alyssa Carlson, and Chris Abraham. Kytopen, Cambridge, MA, USA.

INTRODUCTION

Efficient, simultaneous delivery of multiple classes of nucleic acids presents unique challenges for engineered cell therapies. The Flowfect® non-viral cell engineering technology offers increased payload capacities required for delivering multiple transgenes, limited in the viral context, but critical to advancing the advanced cell therapies. Flowfect™ platforms offer a seamless, scalable continuous flow transfection process, from small volumes (50-100 µL), geared towards discovery and optimization experiments on the Flowfect Discover™ 96-well platform to manufacturing volumes (1 mL- >L) for processing up to billions of cells on the Flowfect Tx™ platform.



Here, we highlight the utility of 1) small-scale, rapid and systematic optimization for efficient co-delivery of mRNA and plasmid DNA and 2) direct scale-up of co-delivery with the Flowfect Tx™ platform in primary human T cells. Initial experiments using the Flowfect Discover™ 96-well platform identified optimal parameters for co-delivery of mCherry mRNA and GFP pDNA in primary activated T cells. These conditions were transferred to the large-scale Flowfect Tx™ platform, achieving comparable results without the need for additional optimization.

To showcase the performance of Flowfect™ platforms in an immunotherapeutic context, we partnered with Bio-technne, a leading provider of cell and gene therapy reagents. Using Bio-technne's TcBuster™ non-viral transposon system, we demonstrate exceptional efficiencies in transposase mediated insertion of a large 5.1kb multi-gene TcBuster™ Transposon CD19CAR-DHFR-eGFP into primary T cells on both Flowfect™ platforms.


KEY TAKEAWAYS

- Leveraging the Flowfect Discover™ 96-well platform, rapid identification and optimization of parameters driving efficient and simultaneous co-delivery (>60%) of mCherry mRNA and GFP pDNA in primary activated T cells were achieved.
- Optimal parameters identified on the Flowfect Discover™ 96-well platform seamlessly transitioned to large-scale volumes on the Flowfect Tx™ GMP platform with no further optimization, achieving high simultaneous co-transfection efficiencies greater than 60% for both mRNA and pDNA while maintaining high viability.
- Supporting the immunotherapeutic relevance, we partnered with Bio-technne, and merging their TcBuster™ transposon system with the Flowfect Discover™ platform achieved up 50% transposition efficiency of a large 5.1 kb multi-gene TcBuster™ Transposon CD19CAR-DHFR-eGFP with >85% viability in primary activated T cells.
- Scaling towards therapeutic requirements, using the Flowfect GMP Tx™ platform, again with Bio-technne's TcBuster™ Transposon CD19CAR-DHFR-eGFP, we achieved up to 55% transposition and >70% cell viability.


TECHNOLOGY



Flowfect Discover™ (Small Volume) Platform Enables Efficient Engineering of CAR-T Cells Using TcBuster™ Non-viral Transposon System



Flowfect Tx™ (Large Volume) Platform Enables Seamless Scale Up for Engineering CAR-T Cells Using TcBuster™ Non-viral Transposon System



CONCLUSIONS

- Harnessing the Flowfect® continuous flow transfection technology drives efficient, simultaneous, co-delivery of multiple nucleic acids. Integrating Bio-technne's TcBuster™ non-viral transposon system, we were able to showcase highly efficient CAR-T engineering while maintaining exceptionally high-viabilities.
- Kytopen's Flowfect® continuous flow transfection technology provides best-in-class performances for both small-scale research utility as well as large-scale manufacturing capability.
- The seamless transition between Kytopen's two Flowfect™ platforms eliminates traditional scale-up bottlenecks while maintaining critical performance attributes, providing a foundation for streamlining cell therapy manufacturing from early development through clinical production.

